



New Knife Mill GRINDOMIX GM 200 for Quick and Effective Sample Homogenization

Tough? Fatty? Fibrous? No problem for the new GM 200! The GRINDOMIX is the perfect mill to homogenize dry, oily, fatty, soft and tough sample materials - for analysis results with minimum standard deviation.

The model GM 200 for sample volumes up to 700 ml has been completely revised. Thanks to the powerful 1000 W drive and the innovative Boost function, the mill is powerful enough to homogenize even difficult samples like tough meat with skin or fibrous plants very quickly and efficiently. The new GM 200 is operated via a user-friendly 4.3" touch display and allows storage of 8 SOPs and 4 program sequences.

The GM 200 is the first RETSCH mill to be controlled via the new RETSCH App. This app includes functions like remote control, display of process parameters or creation and storage of application information.

New features of the GM 200:

- Convenient operation via 4.3" touch display
- Homogenization even of difficult samples within seconds thanks to powerful 1000 W drive
- Boost function with 14.000 rpm facilitates size reduction of fibrous and sticky samples
- Memory for 8 Standard Operating Procedures and 4 program sequences
- Control via optional RETSCH App

Press Release

RETSCH is the leading solution provider for neutral-to-analysis sample preparation and characterization of solids. On the basis of more than 100 years of experience RETSCH develops innovative size reduction and sieving equipment which is characterized by excellent performance, operating convenience, safety and a long lifetime.

www.retsch.com/gm200