

Riplate® medio – more volume, less space!

Due to the compact design Riplate® medio enables more working volume with half the stacking height.

Increased wells to avoid cross-contamination

The tangentially contacting wells of the Riplate® medio leads to more volume. The higher well rims provides a good bioturbation at reduced risk of cross-contamination and allow an airtight closure or sealing with Ritter Capmats.

Like all Ritter Riplates® they offer a high resistance towards commonly used chemicals, solvents and alcohols and also a high stability during centrifugation. Riplate® medio plates are available in sizes 0.5 ml or 1 ml per well and also as sterile version if required.

User benefits summarized:

- Space efficiency with up to 50% lower stack height
- More capacity per well
- Cylindrical wells enable a good bioturbation of samples
- Raised rims reduce risk of cross-contamination
- Airtight closure or sealing
- Alphanumeric coding for data tracking

Two sizes available

The Riplate® medio 0.5 ml plate is 20.5 mm high, offering total volume of 0.6 ml per well with working volume from 20 µl to 500 µl.

The Riplate® medio 1 ml version is 27 mm high, offering total volume of 1.2 ml per well with working volume from 50 µl to 1000 µl.

↖

Artikel-Nr. <i>Article no.</i>	Artikelbezeichnung <i>Article name</i>	VE <i>PU</i>	Stück/Umkarton <i>Pieces/cover. box</i>
43001-0101	Riplate® medio, PP, 1 ml	32	320
43001-0102	Riplate® medio, PP, 0,5 ml	40	400
43001-1101	Riplate® medio, PP, 1 ml, sterile	10	100
43001-1102	Riplate® medio, PP, 0.5 ml, sterile	10	100

Zubehör

Artikel-Nr. <i>Article no.</i>	Artikelbezeichnung <i>Article name</i>	VE <i>PU</i>	Stück/Umkarton <i>Pieces/cover. box</i>
43001-0104	Capmat Riplate® medio	20	100
43001-1104	Capmat Riplate® medio, sterile	10	100

Ritter GmbH

The Ritter GmbH was founded 1965. The company is established in Schwabmünchen, Bavaria. On 25,000 m² Ritter develops and produces with about 300 employees high-quality plastic products for the global market, certified to ISO 9001 and EN ISO 13485 standard.

From product idea to development on modern CAD systems, the in-house tool and mold-production, quality assurance, logistics and distribution is all performed under one roof at Ritter. In more than 50 countries millions of products of Ritter Medical are used in clinical and biotechnological applications every day.



Riplate® medio 1ml and 0,5 ml available



Capmat medio – for an airtight closure or sealing and suitable for both sizes of the Riplate® medio

Source: Ritter GmbH, Schwabmünchen